

Adenine

Catalog No: #S0064



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenine
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	135.13
Formulation	C5H5N5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Metabolism system
SMILES: NC1=NC=NC2=C1NC=N2
Pathways: Neuroscience
Receptor: Adenine phosphoribosyltransferase
Boiling pt: >360°C
Solubility: Ethanol:

References

1. Barrett C, et al. Assay Drug Dev Technol. 2006 Dec;4(6):661-9.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.