

Oxiniacic Acid

Catalog No: #S0067

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxiniacic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	139.11
Formulation	C6H5NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Metabolism systemSMILES:OC(=O)C1=C[N+](O)=CC=C1Pathways:OthersReceptor:OthersBoiling pt:469CMelting pt:254-255CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Atac A, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2012 Jan;85(1):145-54.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.