

Uracil

Catalog No: #S0069

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Uracil
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	112.09
Formulation	C4H4N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OtherSMILES:O=C1NC=CC(=O)N1Pathways:OthersReceptor:RNABoling pt:Melting pt:>300CSolubility:Ethanol: 22 mg/mL (196.27 mM); DMSO: 22 mg/mL (196.27 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Imming P, et al. Nat Rev Drug Discov. 2006 Oct; 5(10):821-34.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.