

Homosalate

Catalog No: #S0071

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Homosalate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	262.35
Formulation	C16H22O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Other
SMILES: CC1CC(CC(C)(C)C1)OC(=O)C1=CC=CC=C1O
Pathways: Others
Receptor: Others
Boiling pt: 341.1°C
Melting pt: Solubility: DMSO: 55mg/mL
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Kim TH, et al. J Toxicol Environ Health A. 2014;77(4):202-13.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.