

Xylitol

Catalog No: #S0589

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Xylitol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	152.15
Formulation	C5H12O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Metabolism systemSMILES:C([C@H](C([C@H](CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 50 mg/mLAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Azarpazhooh A, et al. Cochrane Database Syst Rev. 2011 Nov 9;(11):CD007095.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.