

JW55

Catalog No: #S1807

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JW55
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR inhibitor
Calculated MW	434.48
Formulation	C25H26N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research

Area: Other SMILES: COC1=CC=C(C=C1)C1(CNC(=O)C2=CC=C(NC(=O)C3=CC=CO3)C=C2)CCOCC1 Pathways: Metabolism Receptor: TNKS1; TNKS2
Boiling pt: Melting pt: Solubility: DMSO, Ethanol: 100 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Waaler J, et al. Y Res. 2012 Jun 1;72(11):2822-32.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.