

JSH-23

Catalog No: #S1930

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JSH-23
Brief Description	Inhibitors
Purification	98.00%
Target Name	NF-κB inhibitor
Calculated MW	240.34
Formulation	C16H20N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Immunology
SMILES: CC1=CC(N)=C(NCCCC2=CC=CC=C2)C=C1
Pathways: NF-κB Receptor: NF-κB
Boiling pt: Melting pt:
Solubility: Ethanol: 20 mg/mL (83.21 mM); DMSO: 48 mg/mL (199.71 mM)
Appearance: Gray Spongiform Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kasparikova J, et al. FEBS J. 2014, doi: 10.1111/febs.12711.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.