



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	YL-109
Brief Description	Inhibitors
Purification	99.77%
Target Name	AhR activator
Calculated MW	257.31
Formulation	C14H11NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Other
SMILES: COC1=C/C(=C/2\NC3=CC=CC=C3S2)/C=CC1=O
Pathways: Metabolism
Receptor: AhR
Boiling pt: Melting
pt: Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hiyoshi H, et al. Sci Rep. 2014 Nov 18;4:7095.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.