

JZL195

Catalog No: #S2338

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JZL195
Brief Description	Inhibitors
Purification	98.00%
Target Name	FAAH inhibitor; Lipase inhibitor
Calculated MW	433.46
Formulation	C24H23N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research

Area: Other SMILES: [O-][N+](=O)C1=CC=C(OC(=O)N2CCN(CC3=CC(OC4=CC=CC=C4)=CC=C3)CC2)C=C1 Pathways: Metabolism Receptor: FAAH; MAGL
Boiling pt: Melting pt: Solubility: DMSO: 50 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Long JZ, et al. Proc Natl Acad Sci U S A. 2009 Dec 1;106(48):20270-5.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.