

Description

Product Name	YM-90709
Brief Description	Inhibitors
Purification	98.00%
Target Name	IL Receptor antagonist
Calculated MW	359.42
Formulation	C22H21N3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research

Area:OthersSMILES:CC1(CC2=CC(=C(C=C2C3=CC4=NC5=CC=CC=C5N=C4N31)OC)OC)CPathways:Immunology/InflammationReceptor:IL-5Bolin
g pt:Melting pt:Solubility:Soluble in DMSOAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in
the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Int Immunopharmacol. 2002 Nov;2(12):1693-702.?Morokata T, et al.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.