

JWH-133

Catalog No: #S3584

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JWH-133
Brief Description	Inhibitors
Purification	98.00%
Target Name	Cannabinoid Receptor agonist
Calculated MW	312.49
Formulation	C22H32O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Others SMILES: CCCC(C)(C)C1=CC2=C(C=C1)[C@@H]3CC(=CC[C@H]3C(O2)(C)C)C Pathways: GPCR/G

Protein Receptor: CB2 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bellini G, et al. PKCβⅡ-mediated cross-talk of TRPV1/CB2 modulates the glucocorticoid-induced osteoclast overactivity. Pharmacol Res. 2017 Jan; 115:267-274.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.