

JNK-IN-7

Catalog No: #S3598

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JNK-IN-7
Brief Description	Inhibitors
Purification	99.53%
Target Name	JNK inhibitor
Calculated MW	493.56
Formulation	C28H27N7O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research

Area: Others SMILES: CN(C)C/C=C/C(=O)NC1=CC=CC(=C1)C(=O)NC2=CC=C(C=C2)NC3=NC=CC(=N3)C4=CN=CC=C4 Pathways: MAPK Receptor: JNK1; JNK2; JNK3
Boiling pt: Melting pt: Solubility: DMSO: >24.7mg/ML Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhang T, et al. Chem Biol, 2012, 19(1): 140-154.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.