

## Description

Product Name	YL0919
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor agonist
Calculated MW	334.84
Formulation	C18H23CIN2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:OthersSMILES:Cl.O=C3C=CC=CN3CC2(O)CCN(Cc1ccccc1)CC2Pathways:GPCR/G Protein; NeuroscienceReceptor:5-HT; 5-HT1ABoling pt:Melting pt:Solubility:Soluble in DMSOAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Chen HX, et al. PLoS One. 2013 Dec 18;8(12):e83271.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.