

YYA-021

Catalog No: #S6256



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	YYA-021
Brief Description	Inhibitors
Purification	98.00%
Target Name	gp120/CD4 inhibitor
Calculated MW	317.43
Formulation	C18H27N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Other
SMILES: CC1=CC=C(C=C1)NC(=O)C(=O)NC2CC(NC(C2)(C)C)(C)C
Pathways: Immunology/Inflammation
Receptor: Gp120
Boiling pt: Melting pt: Solubility: DMSO: 0.01 mg/mL (0.03 mM) Ethanol:

References

1. Otsuki H, et al. J Gen Virol. 2013, 94(Pt 12), 2710-2716.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.