

# ZM241385

Catalog No: #S7021



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZM241385
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adenosine Receptor antagonist
Calculated MW	337.34
Formulation	C16H15N7O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Product Description

Research Area: Others  
SMILES: C1=COC(=C1)C2=NN3C(=NC(=NC3=N2)NCCC4=CC=C(C=C4)O  
NPathways: GPCR/G Protein Receptor: A2A  
Boiling pt: Melting pt: Solubility: Ethanol

## References

1. Keddie JR., et al. In vivo characterisation of ZM 241385, a selective adenosine A2A receptor antagonist. Eur J Pharmacol. 1996 Apr 22;301(1-3):107-13.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.