

# Mouse Defensin-6(DEFA6) ELISA kit

Catalog No: #EK20241



Package Size: #EK20241-1 48T #EK20241-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Defensin-6(DEFA6) ELISA kit
Brief Description	ELISA kit
Applications	ELISA
Species Reactivity	Ms

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.126ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.