

Mouse Acetoacetyl-CoA synthetase (AACS) ELISA Kit

Catalog No: #EK20260



Package Size: #EK20260-1 48T #EK20260-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Acetoacetyl-CoA synthetase (AACS) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.9pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.