

Human IgG Fc fragment

Catalog No: #AP89633



Package Size: #AP89633-1 1mg #AP89633-2 10mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human IgG Fc fragment
Brief Description	Human IgG Fc fragment
Other Names	Fragment crystallizable region
Target Sequence	Full length protein
Formulation	PBS, pH 7.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>95% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.