

ZNF364 Rabbit mAb

Catalog No: #52012

Package Size: #52012-1 50ul #52012-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF364 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-1B3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ZNF364
Conjugates	Unconjugated
Modification	Unmodification
Other Names	E3 ubiquitin-protein ligase RNF115; Rabring 7; Zinc finger protein 364
Accession No.	Swiss-Prot:Q9Y4L5GeneID:27246
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of ZNF364/RFN115 in Jurkat, HeLa cell lysates using ZNF364/RFN115 Rabbit mAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.

Background

E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes. Through their polyubiquitination, may play a role in the endosomal trafficking and degradation of membrane receptors including EGFR, FLT3, MET and CXCR4.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.