

EWSR1 Rabbit mAb

Catalog No: #52273

Package Size: #52273-1 50ul #52273-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EWSR1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S04-4B3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human EWSR1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	EWS; EWS-FLI1; bK984G1.4
Accession No.	Swiss-Prot:Q01844GeneID:2130
Calculated MW	Calculated MW: 69 kDa; Observed MW: 85 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000;

Images



Western blot detection of EWSR1 in HeLa cell lysates using EWSR1 Rabbit mAb(1:1000 diluted). Predicted band size:69kDa.Observed band size:85kDa.

Background

Swiss-Prot Acc.Q01844.Might normally function as a transcriptional repressor. EWS-fusion-proteins (EFPS) may play a role in the tumorigenic process. They may disturb gene expression by mimicking, or interfering with the normal function of CTD-POLII within the transcription initiation complex. They may also contribute to an aberrant activation of the fusion protein target genes.SitesFeature keyPosition(s)Description

Note: This product is for in vitro research use only and is not intended for use in humans or animals.