

ORP1 Rabbit mAb

Catalog No: #52479

Package Size: #52479-1 50ul #52479-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ORP1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-2E2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ORP1
Conjugates	Unconjugated
Modification	Unmodification
Accession No.	Swiss-Prot:Q9BXW6GeneID:114876
Calculated MW	Calculated MW: 109 kDa; Observed MW: 109,50 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of ORP1 in Rat Brain,C6,HeLa cell lysates using ORP1 Rabbit mAb(1:1000 diluted).Predicted band size:109kDa.Observed band size: 109,50kDa.

Background

Swiss-Prot Acc.Q9BXW6.Oxysterol-binding protein-related protein 1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.