

TPP1 Rabbit mAb

Catalog No: #52626



Package Size: #52626-1 50ul #52626-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TPP1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S09-2G9
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human TPP1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	CLN2; GIG1; LPIC; SCAR7; TPP-1
Accession No.	Swiss-Prot:O14773GeneID:1200
Calculated MW	Calculated MW: 61,35 kDa; Observed MW: 61,48 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/2000;

Images



Western blot detection of TPP1 in HeLa, Jurkat, HepG2 cell lysates using TPP1 Rabbit mAb (1:1000 diluted). Predicted band size: 61,35 kDa. Observed band size: 61,48 kDa.

Background

Swiss-Prot Acc.O14773. Lysosomal serine protease with tripeptidyl-peptidase I activity. May act as a non-specific lysosomal peptidase which generates tripeptides from the breakdown products produced by lysosomal proteinases. Requires substrates with an unsubstituted N-terminus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.