

## Vitamin D Binding Protein Rabbit mAb

Catalog No: #52648



Package Size: #52648-1 50ul #52648-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Vitamin D Binding Protein Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-5G8
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human Vitamin D Binding protein
Conjugates	Unconjugated
Modification	Unmodification
Other Names	DBP; VDB; GRD3; VDBG; VDBP; GcMAF; DBP/GC; Gc-MAF; DBP-maf; HEL-S-51
Accession No.	Swiss-Prot:P02774GenelD:2638
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/1000-1/5000;

## Images



Western blot detection of Vitamin D Binding protein in Rat Brain, 3T3, HeLa cell lysates using Vitamin D Binding protein Rabbit mAb (1:1000 diluted). Predicted band size: 53kDa. Observed band size: 53kDa.

## Background

Swiss-Prot Acc.P02774. Involved in vitamin D transport and storage, scavenging of extracellular G-actin, enhancement of the chemotactic activity of C5 alpha for neutrophils in inflammation and macrophage activation.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.