

NUP62CL Rabbit mAb

Catalog No: #52682

Package Size: #52682-1 50ul #52682-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUP62CL Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-1A5
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human NUP62CL
Conjugates	Unconjugated
Modification	Unmodification
Other Names	NUP62L; Nucleoporin-62 C-terminal-like protein
Accession No.	Swiss-Prot:Q9H1M0GeneID:54830
Calculated MW	Calculated MW:21 kDa, Observed MW:21 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of NUP62CL in HeLa, Ramos cell lysates using NUP62CL Rabbit mAb(1:1000 diluted). Predicted band size:21kDa. Observed band size:21kDa.

Background

This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore complexes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.