

IKZF3 Rabbit mAb

Catalog No: #52781



Package Size: #52781-1 50ul #52781-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IKZF3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-0C5
Isotype	IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human IKZF3
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AIO; AIOLOS; ZNFN1A3
Accession No.	Swiss-Prot:Q9UKT9GeneID:22806
Calculated MW	Calculated MW:58 kDa,Observed MW:58 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/1000 IHC: 1/20

Images



Western blot detection of IKZF3 in HL-60,U251 cell lysates using IKZF3 Rabbit mAb(1:1000 diluted).Predicted band size:58kDa.Observed band size:58kDa.



Immunohistochemistry of IKZF3 in paraffin-embedded Human tonsil using IKZF3 Rabbit mAb at dilution 1/50

Background

Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.