

NUP155 Rabbit mAb

Catalog No: #52825

Package Size: #52825-1 50ul #52825-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUP155 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S03-6F7
Isotype	IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human NUP155
Conjugates	Unconjugated
Modification	Unmodification
Other Names	N155; ATFB15
Accession No.	Swiss-Prot:O75694GeneID:9631
Calculated MW	Calculated MW:155 kDa,Observed MW:155 kDa
Formulation	50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of NUP155 in HeLa,A549,HL-60 cell lysates using NUP155 Rabbit mAb(1:1000 diluted).Predicted band size:155kDa.Observed band size:155kDa.

Background

Essential component of nuclear pore complex. Could be essential for embryogenesis. Nucleoporins may be involved both in binding and translocating proteins during nucleocytoplasmic transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.