

Mouse Peptidase Inhibitor 16 (PI16) ELISA Kit

Catalog No: #EK20247



Package Size: #EK20247-1 48T #EK20247-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Peptidase Inhibitor 16 (PI16) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Accession No.	Swiss-Prot#:Q9ET66

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.118ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.