

Rat N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit

Catalog No: #EK20248

Package Size: #EK20248-1 48T #EK20248-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.128ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.