

# Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit

Catalog No: #EK20254



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20254-1 48T #EK20254-2 96T

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) ELISA Kit |
| Brief Description  | ELISA Kit  |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Accession No.      | Swiss Prot:Q9CPV9  |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.