

N4BP1 Rabbit Polyclonal Antibody

Catalog No: #53795



Package Size: #53795-1 50ul #53795-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	N4BP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse N4BP1 (NP_085040.2).
Other Names	N4BP1
Accession No.	Swiss Prot:O75113Gene ID:80750
Calculated MW	99kDa
SDS-PAGE MW	120kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using N4BP1 antibody.



Immunofluorescence analysis of U2OS cells using N4BP1 Rabbit pAb.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.