

MCAM Rabbit Polyclonal Antibody

Catalog No: #54322

Package Size: #54322-1 50ul #54322-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MCAM Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human MCAM (NP_006491.2).
Other Names	MCAM;CD146;MUC18
Accession No.	Swiss Prot:P43121Gene ID:4162
Calculated MW	57kDa/71kDa
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded human colon carcinoma using MCAM antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.