

# DNAJC17 Rabbit Polyclonal Antibody

Catalog No: #54644



Package Size: #54644-1 50ul #54644-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DNAJC17 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human DNAJC17.
Other Names	—
Accession No.	Swiss Prot:Q9NVM6Gene ID:55192
Calculated MW	35kDa
SDS-PAGE MW	35KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:2000 - 1:5000 IHC □ 1:50 - 1:200


## Images




Western blot analysis of extracts of various cell lines, using DNAJC17 antibody.



Western blot analysis of extracts of various cell lines, using DNAJC17 antibody.



Immunohistochemistry of paraffin-embedded human esophageal cancer using DNAJC17 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse spinal cord using DNAJC17 Rabbit pAb.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.