

Recombinant Human Deformed epidermal autoregulatory factor 1 homolog(DEAF1)

Catalog No: #AP71005

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP71005-1 20ug #AP71005-2 100ug #AP71005-3 1mg

Description

Product Name	Recombinant Human Deformed epidermal autoregulatory factor 1 homolog(DEAF1)
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	Deaf1; DEAF1 transcription factor; DEAF1_HUMAN; Deformed epidermal autoregulatory factor 1 homolog; Nuclear DEA1 related transcriptional regulator; Nuclear DEAF-1-related transcriptional regulator; NUDR; SPN; Suppressin; Zinc finger MYND domain containing protein 5; Zinc finger MYND domain-containing protein 5; ZMYND5
Tag Info	N-terminal His-tagged
Target Sequence	1-565 Full length protein MEDSDSAAKQ LGLAEAAVA AAAAVAAAAA AAAGGEAEEP VLSRDEESEE DADSEAERET PRVTAVAVMA AEPGHMDMGA EALPGPDEAA AAAFAEVTT VTVANVGAAA DNVFTTSVAN AASISGHVLS GRTALQIGDS LNTEKATLIV VHTDGSIVET TGLKGPAAPL TPGPQSPPTP LAPGQEKGGT KYNWDPSVYD SELPVRCRNI SGTLYKNRLG SGGRGRCIKQ GENWYSPTEF EAMAGRASSK DWKRSIRYAG RPLQCLIQDG ILNPHAASCT CAACDDMTL SGPVRLFPY KRRKKENELP TTPVKKDSPK NITLLPATAA TTFTVTPSGQ ITTSGALTFD RASTVEATAV ISESPAQGDV FAGATVQEAS VQPPCRASHP EPHYPGYQDS CQIAPFPEAA LPTSHPKIVL TSLPALAVPP PTPTKAAPPA LVNGLELSEP RSWLYLEEMV NSLLNTAQQ KTLFEQAKHA STYREAATNQ AKIHADAERK EQSCVNCGRE AMSECTGCHK VNYCSTFCQR KDWKDHQHIC GQSAAVTVQA DEVHVAESVM EKVTV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.