

# Recombinant Human Tumor necrosis factor receptor superfamily member 10D(TNFRSF10D),partial

Catalog No: #AP71091

Package Size: #AP71091-1 20ug #AP71091-2 100ug #AP71091-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Tumor necrosis factor receptor superfamily member 10D(TNFRSF10D),partial
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	CD 264; CD264; DcR 2; DcR2; Decoy receptor 2; Decoy with truncated death domain; OTTHUMP00000123528; TNF receptor related receptor for TRAIL ; TNF related apoptosis inducing ligand receptor 4; TNF-related apoptosis-inducing ligand receptor 4; TNFRSF 10D; TNFRSF10D; TR10D_HUMAN; TRAIL R4; TRAIL receptor 4; TRAIL receptor with a truncated death domain; TRAIL-R4; TRAILR4; TRUND; Tumor necrosis factor receptor superfamily member 10D; Tumor necrosis factor receptor superfamily; member 10d; decoy with truncated death domain
Tag Info	N-terminal His-tagged
Target Sequence	56-211aa Partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

&gt;85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.