

# Recombinant Human Cytochrome P450 4F8(CYP4F8),partial

Catalog No: #AP71137



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP71137-1 20ug #AP71137-2 100ug #AP71137-3 1mg

## Description

Product Name	Recombinant Human Cytochrome P450 4F8(CYP4F8),partial
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	CYP4F8; Cytochrome P450 4F8; CYP4F8
Tag Info	N-terminal His-tagged
Target Sequence	38-520aa Partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.