

Recombinant Mouse Transforming growth factor beta-1(Tgfb1),partial

Catalog No: #AP71114

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP71114-1 20ug #AP71114-2 100ug #AP71114-3 1mg

Description

Product Name	Recombinant Mouse Transforming growth factor beta-1(Tgfb1),partial
Brief Description	Recombinant Protein
Target Name	His-tagged
Accession No.	SwissProt:P04202
Calculated MW	16.9 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	279-390aa Partial ALDTNYCFSSTEKNCCVRQLYIDFRKDLGWKWIHEPKGYHANFCLGPCPYIWSLDTQYSKVLALYNQHNP SASPCCVPQALEPLPIVYYVGRKPKVEQLSNMIVRSCKCS
Formulation	Tris/PBS-based buffer, 5%-50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.