

# Recombinant Mouse Latent-transforming growth factor beta-binding protein 1(Ltbp1) ,partial

Catalog No: #AP71144

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP71144-1 20ug #AP71144-2 100ug #AP71144-3 1mg

## Description

Product Name	Recombinant Mouse Latent-transforming growth factor beta-binding protein 1(Ltbp1) ,partial
Brief Description	Recombinant Protein
Target Name	His-tagged
Other Names	Ltbp1Latent-transforming growth factor beta-binding protein 1; LTBP-1; Transforming growth factor beta-1-binding protein 1; TGF-beta1-BP-1
Accession No.	SwissProt:Q8CG19
Tag Info	N-terminal 6xHis-tagged
Target Sequence	Partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

&gt;85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.