

# General Vitamin C (VC) ELISA Kit

Catalog No: #EK20277



Package Size: #EK20277-1 48T #EK20277-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                  |
|--------------------|----------------------------------|
| Product Name       | General Vitamin C (VC) ELISA Kit |
| Brief Description  | ELISA Kit                        |
| Applications       | ELISA                            |
| Species Reactivity | General                          |

## Application Details

Detect Range:3.12-200 ng/mL

Sensitivity:1.16ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.