

General Homovanillic Acid (HVA) ELISA Kit

Catalog No: #EK20311



Package Size: #EK20311-1 48T #EK20311-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	General Homovanillic Acid (HVA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	General

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.337ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.