

# Gasdermin-D (Full Length+N terminal) Conjugated Antibody

Catalog No: #C53120

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #C53120-AF350 100ul #C53120-AF405 100ul #C53120-AF488 100ul #C53120-AF555 100ul #C53120-AF594 100ul

#C53120-AF647 100ul #C53120-AF680 100ul #C53120-AF750 100ul #C53120-Biotin 100ul #C53120-Conjugated 50ul

## Description

Product Name	Gasdermin-D (Full Length+N terminal) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse GSDMD (NP_081236.1).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GSDMD;DF5L;DFNA5L;FKSG10;GSDMDC1
Accession No.	Swiss Prot:P57764Gene ID:69146
Calculated MW	53kDa
SDS-PAGE MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Background

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.