

EDIL3 Conjugated Antibody

Catalog No: #C54091

Package Size: #C54091-AF350 100ul #C54091-AF405 100ul #C54091-AF488 100ul #C54091-AF555 100ul #C54091-AF594 100ul #C54091-AF647 100ul #C54091-AF680 100ul #C54091-AF750 100ul #C54091-Biotin 100ul #C54091-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EDIL3 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Human, Mouse
Immunogen Description	A synthetic peptide of human EDIL3 (NP_005702.3).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EDIL3; DEL1
Accession No.	Swiss Prot:O43854 GeneID:10085
Calculated MW	52kDa/53kDa
Formulation	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Background

The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.