

Human Kallikrein 1 (KLK1) ELISA Kit

Catalog No: #EK20475



Package Size: #EK20475-1 48T #EK20475-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------------|
| Product Name | Human Kallikrein 1 (KLK1) ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |

Application Details

Detect Range:9.37-600pg/mL

Sensitivity:3.1pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.