

# Recombinant Human Protein max(MAX)

Catalog No: #AP89669

Package Size: #AP89669-1 20ug #AP89669-2 100ug #AP89669-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Protein max(MAX)
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	MAX; BHLHD4Protein max; Class D basic helix-loop-helix protein 4; bHLHd4; Myc-associated factor X
Accession No.	P61244
Target Species	Hu
Tag Info	N-terminal His-tagged
Target Sequence	2-160 Full Length of Mature Protein SDNDDIEVE SDEEQPRFQS AADKRAHHNA LERKRRDHIK DSFHSLRDSV PSLQGEKASR AQILDKATEY IQYMRRKNHT HQQDIDDLKR QNALLEQQVR ALEKARSSAQ LQTNYPSSDN SLYTNAKGST ISAFDGGSDS SSESEPEEPQ SRKCLRMEAS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.