

Recombinant Human Tapasin(TAPBP),partial

Catalog No: #AP89687

Package Size: #AP89687-1 20ug #AP89687-2 100ug #AP89687-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Tapasin(TAPBP),partial
Brief Description	Recombinant Protein
Host Species	Yeast
Target Name	His-tagged
Other Names	TAPBP; NGS17; TAPA; Tapasin; TPN; TPSN; NGS-17; TAP-associated protein; TAP-binding protein
Accession No.	O15533
Target Species	Hu
Tag Info	N-terminal His-tagged
Target Sequence	Partial
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.