

Recombinant Human Serine/threonine-protein kinase LATS1(LATS1) ,partial

Catalog No: #AP89694

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89694-1 20ug #AP89694-2 100ug #AP89694-3 1mg

Description

| | |
|-------------------|--|
| Product Name | Recombinant Human Serine/threonine-protein kinase LATS1(LATS1) ,partial |
| Brief Description | Recombinant Protein |
| Host Species | Yeast |
| Target Name | His-tagged |
| Other Names | h-warts; Large tumor suppressor homolog 1; LATS large tumor suppressor homolog 1; LATS1; LATS1_HUMAN; Serine threonine protein kinase LATS1; Serine/threonine-protein kinase LATS1; WARTS; WARTS protein kinase; wts |
| Accession No. | O95835 |
| Target Species | Hu |
| Tag Info | N-terminal His-tagged |
| Target Sequence | Partial |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.