

Rat IL-1 β ELISA Kit

Catalog No: #EK0532



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat IL-1 β ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range: 62.5-4000 pg/mL

Sensitivity: 15 pg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.