

Recombinant Chicken Indian hedgehog protein(IHH)

Catalog No: #AP89638

Package Size: #AP89638-1 20ug #AP89638-2 100ug #AP89638-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Chicken Indian hedgehog protein(IHH)
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 24-408aa Sequence Info:Full Length of Mature Protein
Target Name	IHH
Other Names	IHH; Indian hedgehog protein; IHH) [Cleaved into: Indian hedgehog protein N-product; Indian hedgehog protein C-product]
Accession No.	Q98938
Target Species	Gallus gallus (Chicken)
Tag Info	N-terminal His-tagged
Target Sequence	CGPGRVV GSRRRPPRKL IPLAYKQFSP NVPEKTLGAS GRYEGKIARN SERFKELTPN YNPDIIFKDE ENTGADRLMT QRCKDRLNSL AISVMNQWPG VKLRVTEGWD EDGHHSEESL HYEGRAVDIT TSDRDRNKYG MLARLAVEAG FDWVYYESKA HIHCSVKSEH SAAAKTGGCF PGRALATLEN GARTPLWALR PGQRVLAMDG AGRPTYSDFL AFLDKEPRAL TAFHVIETRQ PPRRLALTPT HLLFVADNAS APAAQFRPTF ASHVQPGHFV LVAVGSGGLQ PAEVVGVRGR TDVGAYAPLT RHGTLVVDV VASCFAVRE QQLAQMAFWP LRLYHSLGPGVQGDGVHW YSGLLYRLGR MLLPPDSFHP LGAPRAES
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.