

Human Six-transmembrane epithelial antigen of prostate 1 (STEAP1) ELISA Kit

Catalog No: #EK20351



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20351-1 48T #EK20351-2 96T

Description

| | |
|--------------------|---|
| Product Name | Human Six-transmembrane epithelial antigen of prostate 1 (STEAP1) ELISA Kit |
| Brief Description | ELISA Kit |
| Species Reactivity | Human |

Application Details

Detect Range:0.625-40ng/mL

Sensitivity:0.238ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.