

Human CMP-N-acetylneuraminase-beta-galactosamide-alpha- 2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit

Catalog No: #EK20337

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20337-1 48T #EK20337-2 96T

Description

Product Name	Human CMP-N-acetylneuraminase-beta-galactosamide-alpha-2,3-sialyltransferase 1 (ST3GAL1) ELISA Kit
Brief Description	ELISA Kit
Species Reactivity	Human

Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.061ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.