

# Recombinant Human Relaxin receptor 1(RXFP1),partial

Catalog No: #AP89642

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89642-1 20ug #AP89642-2 100ug #AP89642-3 1mg

## Description

Product Name	Recombinant Human Relaxin receptor 1(RXFP1),partial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Sequence Info:partial
Target Name	His-tagged
Other Names	RXFP1; LGR7; Relaxin receptor 1; Leucine-rich repeat-containing G-protein coupled receptor 7; Relaxin family peptide receptor 1
Accession No.	Q9HBX9
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.